FR

DESCRIPTION

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UK

DESCRIPTION

R.T.R. calcium phosphates are manufactured using the latest technologies, which guarantees an extremely high purity level of the products.

R.T.R. cone is made from β tricalcium phosphate granules coated with a matrix of highly purified collagen fibres of bovine origin, which meet the health and safety requirements recommended by the WHO and the EU.

PROPERTIES

R.T.R. is available as granules with a size ranging from 500 μ m to 1 mm.

The size of macropores varies from 100 μ m to 400 μ m and micropores are less than 10 μ m in diameter. These specific structural properties make it possible for newly formed bone to colonise macropores.

R.T.R. is gradually resorbable.

INDICATIONS

R.T.R. cone is designed for the replacement of bone losses, after being impregnated with the patient's blood or with physiological solution, for the filling of tooth sockets for instance.

R.T.R. cone is indicated particularly when the socket cannot be closed, since granules would probably leak out of the socket.

CONTRAINDICATIONS :

- Acute osteomyelitis or necrosis of the recipient site
- Bone degeneration
- Known active infection
- Known allergy to any constituent of the bone substitute

WARNINGS AND PRECAUTIONS FOR USE

The manufacturer of R.T.R. guarantees the material and manufacture quality.

Several criteria are of crucial importance for the operation with R.T.R. to be successful:

- appropriate selection of the patient who will undergo the operation
- rigorous asepsis

The practitioner is responsible for any complication which may result from improper use, faulty operative technique or lack of asepsis. These complications cannot be attributed to the manufacturer.

Underlying oral pathologies such as infections should be treated and ruled out before applying R.T.R..

Avoid saliva contamination.

The bone site which will be in contact with R.T.R. should be debrided and, if possible, well vascularised, so as to create a high quality recipient bed.

Any non-used part of R.T.R. should be systematically disposed of since the product is no longer sterile.

Check the packaging integrity before use.

Do not use if packaging has been damaged or opened before first use.
Do not resterilise.

INSTRUCTIONS FOR USE

R.T.R. cone is inserted into the dental socket after waiting that the socket has been filled with blood. It is always better to suture the operative site. However, should suturing not be possible, make sure that the cone has been completely covered with and protected by the wound edges.

STORAGE

Keep only in the original container.

EXPIRY DATE

Do not use after the expiry date mentioned on the outer package.

PRESENTATIONS

2 x 0,3 cm³ cone (Ø 6mm, H 10mm) in single-unit package sterilised by gamma radiation (minimal dose: 25 kGy).

For professional dental use.

Distributor: SEPTODONT Ltd.

Units R & S - Orchard Business Center
St Barnabas Close Allington, Maidstone
Kent ME16 0JZ - UK

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