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(EN)

Before using a RTR bone substitute, the surgeon must take cognisance of the following:

GENERAL INFORMATION

The manufacturer of RTR bone substitutes guarantees the quality of manufacture and material.

Several criteria are of decisive importance for the success of an intervention using a bone substitute:

- The appropriate selection of the patient for the intervention in question
- The strictest asepsis

The surgeon is responsible for the complications that may result from an erroneous indication and a faulty operating technique or fault involving aseptis, and these cannot be imputed to the manufacturer.

MANUFACTURING - MATERIALS

The RTR bone substitutes are made of tricalcium phosphate β (TCP; Ca3(PO4)2; origin: synthetic) according to the latest technologies to ensure an excellent level of purity. However, their mechanical and functional characteristics can only be guaranteed if the instructions in this manual have been observed. The combination of the chemical composition (pure TCP) and

specific structural properties allows colonisation of macropores by neoformed bone and progressive but total resorption of RTR implants: the macropores have a size that varies from 100 µm to 400 µm for an interconnection rate greater than 80%).

STORAGE AND HANDLING

- Bone substitutes must be stored in their original packaging.
- Bone substitutes should not be modified in any manner other than those authorised by this leaflet.
- Damaged bone substitutes may not be implanted.
- Do not exceed the expiry date indicated on the packaging.

PACKAGING AND STERILISATION

Bone substitutes are delivered sterile.

- They are sterilised with gamma radiation at a minimum dose of 25 kGy.
- The integrity of the individual sterility protector should be checked before use.
- Do not use if the packaging is damaged ☹

RE-STERILISATION AND REUSE

The re-sterilisation and/or reuse of RTR bone substitutes is strictly forbidden, as this may expose the patient to contamination that may lead to infection. The re-sterilisation and/or reuse of a RTR bone substitute would result in a decrease in the performance of the device, namely reduced mechanical strength (risk of fracture) and an incomplete bone remodelling phenomenon. There is also a risk of a significant inflammatory reaction (for example, a host reaction against the graft if it has been implanted before).

ELIMINATION AFTER USE

In the case of extraction of the device, please handle the components in accordance with the centre's medical waste disposal procedure.

INDICATIONS

RTR bone substitutes are intended for filling and reconstructing bone loss.

The choice of the implant is up to the surgeon. Once the manufacturer's recommendation not to place the products under mechanical stress without osteosynthesis material is followed, the indications described below can be addressed using RTR bone substitutes. :

- Bone defects and cavities requiring filling
- Alveolar filling post-extractionally
- Implant reinforcement
- Elevation of the sinus floor

PRESENTATIONS AVAILABLE:

- Curved 0.8 cm³ syringe packaged in an individual package sterilised with gamma radiation (minimum dose 25 kGy).
- Bottle containing 2 cm³ of granules packaged in individual packaging sterilised with gamma radiation (minimum dose of 25 kGy).

RECOMMENDATIONS

The surgeon must choose the size and shape of the graft best suited to the surgical treatment. The radiological examination will allow him to make this choice.

The bone site intended to be in contact with the graft must be abraded and, if possible, well vascularised in order to have a good blood supply.

- In order for the graft to function properly, it must be placed in close and stable contact with the bone site (by compression or osteosynthesis system). In addition, to improve the graft, it is strongly recommended to impregnate the implant with bone marrow or, failing that, autologous blood and to add, where possible, spongy chips from a piercing, a milling or even an abrasion.

Bone substitutes should be used in areas of mechanical stress only if they are used with osteosynthesis equipment that is capable of ensuring the stability of the surgical setup, although this is the case for any reconstruction with bone of autologous or allogeneic origin.

For alveolar filling, the bone substitute is placed in the alveolus taking care to wait until it is filled with blood. Surgical site suturing is always preferable; however, in the case where it is not feasible, care must be taken to recover and protect the bone substitute by the edges of the wound.

Care must be taken to avoid salivary contamination.

RTR bone substitutes are compatible with MRI without artifact risk.

Bone substitutes should be used with caution in patients treated with bisphosphonates and in specific clinical situations (tumor, ongoing chemotherapy, immunodeficiency, poorly controlled diabetes, allergies, irradiated bone).

- The placement should be done without excessive compression of bone substitutes in the cavity to be filled.
- The safety of bone substitutes during pregnancy and lactation has not been established. It is therefore recommended to avoid the use of these bone substitutes during this period, when the condition of the patient allows it.
- Bone substitutes allow several modes of application:
 - In the case of the syringe, the user can apply the granules directly into the site of the bone defect before soft tissue closure. To do this, draw blood or saline through the mesh of the tip of the syringe by pulling on the plunger, until all the granules are soaked. Excess liquid can be removed by pressing the plunger. Remove the tip and place the contents in the surgical site by pressing the plunger. The granules can also be mixed with the patient's blood or saline in a sterile cup, and then placed on the site with a dental spatula.
 - In the case of bottles, the granules can be mixed with the patient's blood or physiological serum in a sterile scoop and then placed on site with a dental spatula. In the case of significant bone defect, the bone substitutes can be mixed with bone particles from the same patient (autologous bone).

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- General contraindications against grafts
- Acute osteomyelitis or necrosis of the recipient site
- Bone degeneration
- Known active infection
- Allergy to a constituent of the bone substitute
 - Where any surgery is contraindicated

COMPLICATIONS

The following complications have sometimes been observed, and these require special attention:

- Delayed healing
- Fracture
- Delayed union and pseudarthrosis
- Pain
- Temporary paralysis of a nerve
- Deep vein thrombosis and pulmonary embolism
- Infection
- Inflammation
- Oedema of the wound
- Allergic reaction
- Opening of the wound
- Reaction to a foreign body
- Tissue degeneration
- Necrosis
- Mobilisation of the dental implant held in place (in whole or in part) by the bone substitute
- Benign nosebleeds (benign epistaxis)
- Benign positional paroxysmal vertigo
- Postoperative seroma
- Failure of bone grafting
- Non-resorption of the bone substitute

EXPCTED PERFORMANCE

Because of a chemical composition that is very close to that of the natural mineral bone, calcium phosphates and in particular tricalcium phosphate β play a significant role in the biological processes governing the stability and regeneration of bone tissue. Thus, the role of RTR bone substitutes is to provide a temporary bone replacement function in the case of significant tissue deficiency of traumatic or physiological origin. These devices thus make it possible to initiate the bone remodelling phenomenon while facilitating their own replacement with newly formed natural bone.

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