

Prof. Joseph Sabbagh (DDS, MSc, PhD, HDR, FICD)

Prof. Sabbagh graduated from Saint Joseph University of Beirut and completed his post-graduate curriculum at the Catholic University of Louvain (UCL), Belgium. Here, he earned his PhD in Biomaterials in 2004 and his Master's in Operative Dentistry (Restorative Dentistry and Endodontics) in 2000. In 2020, he obtained his Higher Degree of Research from the Lebanese University, Beirut.

Prof. Sabbagh is currently a Professor at the Lebanese University, Beirut, where he has directed several research projects and previously acted as Head of the Department of Restorative Dentistry and Endodontics. His private practice focuses on aesthetic dentistry and endodontics in Belgium and Beirut.

Prof. Sabbagh has published more than 25 papers in international peer-reviewed dental journals and has lectured locally and internationally. He has edited and published a book by Springer entitled *Bulk Fill Resin Composites in Dentistry*. He is a member of several dental associations and boards.



Demystifying root canal treatments in daily practice using GenENDO instruments and BioRoot™ Flow*

Introduction

Root canal treatment (RCT) is a common endodontic procedure aiming to preserve a tooth through proper cleaning, shaping and obturation of the root canal system. When properly performed, root canal treatment has an estimated success rate of 90.3% (Ricucci et al., 2011). For retreatments, the success rate is between 65.5% and 77.6% (Stueland H et al., 2023).

Recent innovations in endodontics are revolutionizing root canal therapy by integrating advanced imaging, precision instruments and enhanced disinfection methods, making the treatment easier and more predictable. Three-dimensional imaging—particularly cone-beam computed tomography (CBCT)—enables detailed visualization of complex canal morphologies and



guided access. Cleaning and shaping of a root canal system relies on using flexible nickel-titanium (NiTi) files, simplified sequences, and improved alloys, together with apex locators and enhanced rotary motors, to facilitate efficient shaping of curved canals while minimizing procedural errors. This enables subsequent steps, including coronal flaring, creating a reproducible glide path, determining the working length, and carrying out biomechanical preparation and disinfection of the root canal system. Ultimately, the procedure aims to achieve effective obturation, ensure the longevity of the restoration, and preserve as much of the natural tooth structure as possible (Elmatary et al., 2025 ; ESE, 2006).

The dental community has understood that an adequate endodontic treatment does not mean enlargement of the canals, but rather a proper disinfection of the root canal system combined with a three-dimensional seal. Several techniques using Gutta-Percha, delivered in different modes,

are used for the obturation phase. For several years, the warm vertical technique, described by Schilder in 1972, was considered the gold standard in endodontics. Although showing good clinical outcomes, this technique is complex and involves several steps. Recent development in biomaterials has led to a new category of products called bioceramics, based mostly on calcium silicate (CSCs). They can be used as cements in restorative dentistry for pulp vitality preservation, or as sealers in endodontics (Dong and Xu, 2023).

During the obturation step, the use of bioceramic sealer combined with a single Gutta-Percha cone makes this phase easier and more predictable.

The following clinical case report details a root canal treatment of an upper molar. The root canal preparation was completed with two NiTi files in continuous rotation (Revo-S+, GenENDO, Septodont). The obturation was done using a recent bioceramic sealer (BioRoot™ Flow, Septodont).

Case presentation

Clinical signs and symptoms

A 59-year-old male patient presented to the dental clinic to complete a root canal treatment after severe spontaneous pain localized at the upper right posterior quadrant. The pain had persisted for five days and intensified with hot stimuli. Consequently, the patient went to a hospital offering an emergency dental service. They created an access cavity in tooth #17 to relieve the pain, placed a temporary restoration, and asked him to continue the treatment with an endodontist.

A preoperative periapical radiograph revealed an apical lesion on the palatal root of tooth #17 (upper right second molar) (*Fig. 01*) with a widened periodontal ligament space around the mesial root. No swelling or sinus tract was visible.

Diagnosis

Based on the patient explanations, clinical and radiographic findings, the tooth was diagnosed with a necrosis following an acute irreversible pulpitis. A non-surgical root canal treatment was



Fig. 01 - Preoperative radiography of tooth #17.

planned to complete the emergency procedure initiated at the hospital.

Procedure and treatment

Following administration of local anesthesia (Septanest, 1:200.000, Septodont, Saint-Maur-des-Fossés, France), the temporary cement was

removed using a diamond bur and a carbide tungsten bur was used to finalize the access cavity. Three canal orifices were identified: mesiobuccal (MB), distobuccal (DB) and palatal (P). The working field was isolated using a latex rubber dam fixed with a Softclamp (Kerr, Orange, USA) to avoid any gingival fluid contamination or accidental instrument swallowing.

A manual stainless steel 010 K-File (GenENDO, Septodont) was used to scout the canals' patency (Fig. 02). Final working length was determined using an apex locator and confirmed radiographically.

Cleaning and shaping were initiated using a crown-down technique with rotary files GenENDO Revo-S+ (Septodont) SC2 and SU. There was no need for coronal flaring due to the size and the straight form of the canals.

The shaping of the canal was done with two NiTi Files: GenENDO Revo-S+ SC2 (25/0.4 symmetrical design) was used first, centered downward on the canal in two or three motions. Once the length was obtained with the SC2, a final shaping was done with the GenENDO Revo-S+ SU (25/.06 asymmetric design) for a uniform taper and

optimal preparation (Fig. 04) following one centered downward and one upward movement with selective wall support, finishing at working length.

Copious irrigation with 3.5% sodium hypochlorite was performed throughout instrumentation, and canals were irrigated with EDTA on a final flush flow to remove the smear layer. Rechecking of canal patency was regularly ensured using a GenENDO K-File 010 between each instrument. Figure 5 shows the access cavity with the three cleaned and prepared canals.

In the absence of swelling and after complete drying of the canals, contemporary endodontic strategy recommend the obturation of the canal system during the same session. After final irrigation, the canals were dried using one or two paper points to avoid over-drying. BioRoot™ Flow (Septodont), a recently introduced bioceramic sealer, was injected into each canal with low pressure (Fig. 06).

Using a single-cone technique, one master cone of Gutta-Percha, previously calibrated and confirmed radiographically, was inserted into each of the canals. The Gutta-Percha points were then cut using a heat-cutting device and compacted with



Fig. 02 - Canal exploration using a manual stainless steel file: GenENDO K-File 010.



Fig. 03 - Canal preparation using NiTi rotary file: GenENDO Revo-S+ SC2 (size 25/taper .04).



Fig. 04 - Final shaping using NiTi rotary file: GenENDO Revo-S+ SU (size 25/taper .06).




Fig. 05 - Access cavity showing the three cleaned and prepared canals.




Fig. 06 - Injection of Bioroot™ Flow in the canals.

an endo plugger. Figure 7 shows the entrance of the canals after obturation and compaction. A Teflon pellet was placed in the access cavity and covered by a temporary cement. Figure 8 shows the postoperative radiography of the canal system properly shaped, cleaned and homogenously filled to the radiographic apex.




Fig. 07 - Canal entrances after obturation and compaction. (Gutta-Percha is cut using a heating device, then compacted with an endo plugger.)

Follow-Up and Outcome

At the one-month follow-up, the patient reported no pain or discomfort. Clinical testing showed no tenderness, and radiographic examination revealed resolution of the periodontal ligament widening and signs of periapical healing. The tooth remained functional and symptom-free.




Fig. 08 - Postoperative radiography showing the shaping and obturation of the three canals.

Discussion

Endodontic treatment of molars can be challenging due to variations and complexities in root canal morphology. Use of magnification, electronic apex location, and simplified rotary instrumentation enhanced the precision of canal shaping. The obturation phase was performed with a bioceramic sealer and a single Gutta-Percha cone, which is being used more and more in endodontics.

Successful root canal treatment depends on adequate diagnosis, canal disinfection, and hermetic obturation. Early intervention prevented the progression of periapical disease and preserved the natural tooth structure, avoiding extraction.

Conclusion

This case highlights the importance of comprehensive diagnostic and clinical protocols in the endodontic treatment of maxillary molars. With proper technique

and patient compliance, even complex molars can be predictably treated, resulting in long-term tooth retention and functional restoration.

References

1. Ricucci D, Rutberg M et al, 2011, A prospective cohort study of endodontic treatments of 1,369 root canals: results after 5 years. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 112(6) : 825-42.
2. Stueland H, Østravik D, Handal T. 2023 Treatment outcome of surgical and non-surgical endodontic retreatment of teeth with apical periodontitis. Int Endod J ; 56 (6) : 686-696.
3. Elmatary, A., Moawad, E., Heidarifar, O. et al. Endodontic access cavity preparation: challenges and recent advancements. Br Dent J 238, 469–475 (2025).
4. Dong X & Xu X. 2023 Bioceramics in Endodontics : Updates and future perspectives. Bioengineering: 10. 354.

*Disclaimer

The views, opinions, and statements expressed in this content are solely those of Dr. Joseph Sabbagh.

Dr. Joseph Sabbagh is solely responsible for the scientific and medical positions presented. The pictures and radiographs used in this case study are the propriety of Dr. Joseph Sabbagh.

The author reports a professional and/or financial relationship with Septodont relevant to this article.

The products referenced in this case study are medical devices. BioRoot™ Flow is available in the United States and the European Union and complies with applicable regulatory requirements (FDA 510(k) or MDR 2017/745). GenEndo Range complies with the applicable requirements of the European Union (MDR 2017/745) and is available in EU only.

Please read the instructions for use carefully before use, paying particular attention to warnings, precautions, and contraindications.

These medical devices are regulated health products. Communication intended for healthcare professionals only.

Septodont

58 rue du Pont de Créteil - 94100 Saint-Maur-des-Fossés - France

www.septodont.com

Follow us on social media channels:

